



## Jupiter & Jupiter SW

AUTOCLAVABLE  
FERMENTERS/  
BIOREACTORS  
(JACKETED &  
SINGLEWALL)



**SOLARIS**  
BIOTECH SOLUTIONS

# SINGLE & PARALLEL FERMENTERS/BIOREACTORS

## JUPITER

The **JUPITER** platform offers multiple autoclavable vessel sizes and designs from 2 up to 10 L total volume. Various aspect ratios and thermoregulation designs are also available. The system is highly configurable, built with high quality components, and offered at a competitive price with no strings attached.

Jupiter is available both jacketed and single-wall (**Jupiter SW**).

**JUPITER** typical applications includes the following:

- Education & Basic research
- Scale-up and scale-down studies
- Process development and optimization

**JUPITER** can be used for:

- Biopharmaceutical
- Biofuels
- Food industry
- Bioremediation
- Bioplastic
- Cosmeceutical
- Nutraceutical



**WHY TO  
INVEST  
IN THIS PRODUCT**

The best ratio  
**Quality/  
Capability/Price**  
on the market

**Parallel control**  
up to 24 units

# SINGLE & PARALLEL FERMENTERS/BIOREACTORS

# JUPITER

## Benefits

Up to 24 units managed with one HMI with innovative PARALLEL process control LEONARDO: smart controller designed to provide a high level of automated management of the fermentation/cultivation processes  
Batch, Fed batch or continuous processes

Powerful/ Accurate **brushless motor**, from 1 to 2000 RPM.  
Online absorbed Torques (Nm) and Power (W) measurements obtaining an indirect density indication of the culture broth.

Modbus Digital sensors



LEDA safe sterile sampling system

Safety: pressure relief valve included in each unit

Compact and modular PCS

Additional parameter in modular external boxes for future PCS upgrade including dCO<sub>2</sub>, cell density, weight, peristaltic pumps, ect




N.4 assignable Watson Marlow pumps in entry level

Wide range of options, 5 different volumes and 2 different ratio H/D

Jacketed (fully removable and cleanable) or single wall, with heating blanket and cooling finger (Jupiter SW)



Different gas mixing strategies with up to 5 TMFC



24" touch HMI

Remote access via PC, tablet/smartphone  
Remote control for after sale assistance



# SINGLE & PARALLEL FERMENTERS/BIOREACTORS



## Modbus Digital sensors

### Why a digital sensor?

Digital sensors (including Cell Density products) have been integrated to the Solaris PCS and Leonardo controlling software, giving the user many benefits over traditional analog sensor outputs. Such benefits include a robust communication protocol not susceptible to signal loss, in-software sensor diagnostic information, parallel calibration/batch calibrations and more.



**Sensor life  
traceability**

**Reducing  
background noise**

## GAS MIXING

Hardware and software adaptability are key to enable the best aeration strategy for each process. Thermal mass flow controllers (TMFC) allow precise flow rate control of individual gasses. Up to 5 TMFC's can be configured within each PCS cube and integrated to the controlling software. The powerful software and control platform allows precise cascade adjustment of multiple parameters to manage gas transfer, OTR, kLa, etc.

- n.1 TMFC included in "entry" level system; additional available as optional.
- Various agitator and baffle designs available
- Automatic gas mixing algorithms
- Toro, sintered and other spargers available



# SINGLE & PARALLEL FERMENTERS/BIOREACTORS



## LEDA sterile sampling system

Technical specifications	
Material	VALOX resin (external) silicone (internal)
Autoclavable	121-133°C (up to 30 minutes)
Residual volume	0.04 mL
Flow rate	165 mL/minute



- Sterile single use sampling system up to 180 sterile sampling per batch.
- Needlefree connector is designed to reduce the risk of contamination during sampling.
- The sterile combination of a syringe (3-5-10-30 ml) and a non return valve guarantees the sterility after sampling until the next use

## Smart PCS



Solaris new modular product design strategy decreases time to market and the number of unique parts in the product architecture, increasing the number of product variants. The result is a lean, flexible and smart PCS, which can be stacked in case of parallel processes through a dedicated support.




Additional parameters in modular external boxes for future PCS upgrade including dCO<sub>2</sub>, Cell Density, Weight, Peristaltic pumps, ect.

## Leonardo 3.0

### USER-FRIENDLY SOFTWARE

Solaris controlling software offers a simply laid out, yet powerful platform for experimental design planning and process control. The graphical user interface enables the intuitive selection and adjustment of control functions. Extracted data is compatible with Window Excel but, in addition, Solaris offers a platform where fermentation data can be easily exported in real time and thus managed. This software is included in the supply and can be installed on an unlimited number of the client's PC or laptops.




Parallel synoptic.

### Do it parallel: smarter..faster

Leonardo allows intuitive and time-saving parallel operations. Up to 24 independent fermentations/cultivations can be carried out simultaneously.

### Do it wireless!

Increase mobility: users have the option to access the platform remotely, via PC, tablet, phone. Remote access is multi-level password protected.



# SINGLE & PARALLEL FERMENTERS/BIOREACTORS

# JUPITER

## Data sheet

Vessel					
Solaris Code	Jupiter 2.0	Jupiter 4.0	Jupiter 6.5	Jupiter 8.0	Jupiter 10.0
Production Code	jpt110300	jpt130395	jpt160395	jpt160480	jpt180480
Total Volume (L)	2,00	4,00	6,50	8,00	10,00
Ratio D/H	1:3,0	1:3,0	1:2,5	1:3,0	1:3,0
Min. Working Volume (L)	0,35	0,60	1,10	1,10	1,60
Max. Working Volume (L)	1,40	2,80	4,50	5,50	7,0
Max. temperature	70°C				
Operating pressure	< 0.5 bar				
Headplate Ports (n.10 in Jupiter 2.0; n.13 in the others)	10: n. 1 Agitation Group, n.1 Gas Sparger, n.1 Gas Overlay, n.1 Gas Out/Condenser, n.1 Sampling/Harvesting, n.1 Temperature, n.1 Multifeed, n.2 Sensors DN12, n.1 Spare. 13: n.1 Agitation Group, n.1 Gas Sparger, n.1 Gas Overlay, n.1 Gas Out/Condenser, n.1 Sampling/Harvesting, n.1 Sterile Sampling System, n.1 Temperature, n.1 Multifeed, n.2 Sensors DN12, n.3 Spare.				
Design	Borosilicate Glass Jacketed Vessel				
Materials	Borosilicate Glass and AISI 316 L				
Sensors length (mm)					
pH	325	425	425	425	425
dO <sub>2</sub>	325	425	425	425	425
Dimensions for autoclave (with Condenser)					
Height (mm)	610	705	705	790	790
Diameter (mm)	275	285	315	315	335
Stirring					
Drive	Brushless Motor				
Speed (rpm)	1-1900	1-1800	1-1700	1-1700	1-1700
Nominal Torque (Nm)	0,9	0,9	0,9	1,1	1,1
Impellers	Select from: Rushtons impellers, Marine Impellers, Pitched blade				
Thermoregulation					
Control	PID Control - Accuracy 0,1 °C - Jacketed with n. 2 Electric Cartridge Heaters and cooling valve				
Total Heater Power (W)	400	600	700	700	700
Gas Control & Gas Mixing					
Sparger and overlay Gas Control	TMFC				
Gas Mixing (Air,CO <sub>2</sub> ,O <sub>2</sub> ,N <sub>2</sub> )	n.1 TMFC (included in entry level) + n.4 solenoid valves or + n. of additional TMFC (up to n.4)				
Sparger type	Select from: Toro type (ring), sintered microbubbling - both provided with 0,22 µm sintered filter				
Gas Out	n. 1 Condenser + 0,22 µm sinterized filter				
Peristaltic Pumps					
	n.4 Watson Marlow type 114, fixed speed, max. 60 rpm, volumetric flow 0,5-51 ml/min, function assignable from software				
	(optional) Watson Marlow type 313 FDM/D, max. speed 350 rpm, volumetric flow 1,5-1750 ml/min, function assignable from software				
Controller					
Master Control Module	From 1 to 24 units - 35x37x36 cm				
HMI with Leonardo software	Operate interface 58x15x48 cm with 24" monitor				

## Controls

INTEGRATED IN THE PCS	Temperature		
	Sensor	PT100	
	Accuracy	0,1 °C	
	Control system	Measuring resident in Leonardo 3.0 software	
	Control range	0 - 70°C	
	pH	Sensor	Digital sensor
		Sensitivity	57 to 59 mV/pH
		Control system	Measuring resident in Leonardo 3.0 software
		Control range	0 - 14
		Operation temperature	0 - 130°C
	Pressure range	0 - 6 bar	
	dO <sub>2</sub>	Sensor	Digital Optical sensor
Accuracy		±0.05%-vol, 21±0.2%-vol, 50±0.5%-vol	
Control system		Measuring resident in Leonardo 3.0 software	
Control range		0,05 - 300% air saturation	
Operation temperature		-10 - 130°C	
Pressure range	0 - 12 bar		
Antifoam/Level	Sensor	Solaris sensor	
	Control	Measuring resident in Leonardo 3.0 software	
Redox (ORP)	Sensor	Digital sensor	
	Sensitivity	57 to 59 mV/pH	
	Control system	Measuring resident in Leonardo 3.0 software	
	Control range	±2000 mV	
	Operation temperature	- 10 -130°C	
Pressure range	≤ 6 bar		
Conductivity	Sensor	Digital sensor	
	Accuracy	±3%	
	Control system	Measuring resident in Leonardo 3.0 software	
	Control range	1 - 3000 µS/cm	
	Operation temperature	0 -130°C	
Pressure range	0 - 20 bar		
dCO <sub>2</sub>	Sensor	Analog sensor	
	Accuracy	±10% (pCO <sub>2</sub> 10-900 mbar) ≥ ±10%(pCO <sub>2</sub> > 900 mbar))	
	Control system	Measuring resident in Leonardo 3.0 software	
	Control range	0,00-200% saturation	
	Operation temperature	-20,0-150°C	
Cell density	Sensor	Digital sensor	
	Accuracy	Mammalian cells in suspension ±5·10 <sup>4</sup> cells/ml - Fermentation ±0.05 g/l dry weight	
	Control system	Measuring resident in Leonardo 3.0 software	
	Pressure range	0-3 bar (option 1) 0-10 bar (option 2)	
	Operation temperature	0-60°C (option 1) 0-80°C (option 2) (max. sterilization temperature 135°C)	
Option 1	Dencytee: Total cell density based on turbidity (Two ranges: 10 <sup>4</sup> 5 to 10 <sup>8</sup> mammalian cells/ml - 0.5 to 100 g/L dry weight)		
	Incyte: Viable cell density based on capacitance (Two ranges: 5x10 <sup>4</sup> 5 to 8x10 <sup>8</sup> mammalian cells/ml - 5 to 200 g/L dry weight)		
EXTERNAL MODULAR BOX	Weight		
	Sensor	Digital balance	
	Accuracy	±0.2 g	
Peristaltic pumps	Control	Measuring resident in Leonardo 3.0 software	
	WM 114	10-60 rpm	

## Chiller

- Optionally JUPITER can be equipped with a chiller for heat removal from your culture minimizing lab water usage
- Using this system you don't need a water supply line in your lab
- Cost-effective cooling of fermenters
- Easy operation
- Refregerant level monitoring



Chiller data sheet	
Working temperature range	-10°C / +40°C
Temperature stability	±0.5
Power consumption	0.7 kW
Filling volume range	2-8 L
Cooling output at 20°C measured with ethanol	0.25-0.60 kW
Cooling output at 10°C measured with ethanol	0.20-0.50 kW
Cooling output at 0°C measured with ethanol	0.15-0.36 kW
Cooling output at -10°C measured with ethanol	0.09-0.15 kW
Pump pressure max.	0.35-1.30 bar
Pump flow max.	16-35 L/min.


# SINGLE & PARALLEL FERMENTERS/BIOREACTORS

# JUPITER SW



## Data sheet

Vessel					
Solaris Code	Jupiter SW 2.0	Jupiter SW 4.0	Jupiter SW 6.5	Jupiter SW 8.0	Jupiter SW 10.0
Production Code	L110300	L130395	L160395	L160480	L180480
Total Volume (L)	2,00	4,00	6,50	8,00	10,00
Ratio D/H	1:3,0	1:3,25	1:2,50	1:3,20	1:3,0
Min. Working Volume (L)	0,35	0,60	1,10	1,10	1,60
Max. Working Volume (L)	1,40	2,80	4,50	5,50	7,0
Max. temperature	70°C				
Operating pressure	< 0.5 bar				
Headplate Ports (n.10 in Jupiter 2.0; n.13 in the others)	10: n. 1 Agitation Group, n.1 Gas Sparger, n.1 Gas Overlay, n.1 Gas Out/Condenser, n.1 Sampling/Harvesting, n.1 Temperature, n.1 Multifeed, n.2 Sensors DN12, n.1 Cooling Finger. 13: n.1 Agitation Group, n.1 Gas Sparger, n.1 Gas Overlay, n.1 Gas Out/Condenser, n.1 Sampling/Harvesting, n.1 Sterile Sampling System, n.1 Temperature, n.1 Multifeed, n.2 Sensors DN12, n.1 Cooling Finger, n.2 Spare.				
Design	Borosilicate Glass Vessel				
Materials	Borosilicate Glass and AISI 316 L				
Sensors length (mm)					
pH	325	425	425	425	425
dO <sub>2</sub>	325	425	425	425	425
Dimensions for autoclave (with Condenser)					
Height (mm)	610	705	705	790	790
Diameter (mm)	275	285	315	315	335
Stirring					
Drive	Brushless Motor				
Speed (rpm)	1-1900	1-1800	1-1700	1-1700	1-1700
Nominal Torque (Nm)	0,9	0,9	0,9	1,1	1,1
Impellers	Select from: Rushtons impellers, Marine Impellers, Pitched blade				
Thermoregulation					
Control	PID Control - Accuracy 0,1 °C - n. 1 Electric Heating Blanket, n.1 cooling finger				
Total Heater Power (W)	100	125	125	160	180
Gas Control & Gas Mixing					
Sparger and overlay Gas Control	TMFC				
Gas Mixing (Air,CO <sub>2</sub> ,O <sub>2</sub> ,N <sub>2</sub> )	n.1 TMFC (included in entry level) + n.4 solenoid valves or + n. of additional TMFC (up to n.4)				
Sparger type	Select from: Toro type (ring), sintered microbubbling - both provided with 0,22 µm sintered filter				
Gas Out	n. 1 Condenser + 0,22 µm sinterized filter				
Peristaltic Pumps					
	n.4 Watson Marlow type 114, fixed speed, max. 60 rpm, volumetric flow 0,5-51 ml/min, function assignable from software				
	(optional) Watson Marlow type 313 FDM/D, max. speed 350 rpm, volumetric flow 1,5-1750 ml/min, function assignable from software				
Controller					
Master Control Module	From 1 to 24 units - 35x37xh36 cm				
HMI with Leonardo software	Operate interface 58x15xh48 cm with 24" monitor				



# SOLARIS

## BIOTECH SOLUTIONS

**SOLARIS BIOTECHNOLOGY srl**

Via Bachelet, 58 - 46047 Porto Mantovano  
Mantova - Italy  
Phone: +39 0376 408760  
Fax: +39 0376 385108  
Email: info@solarisbiotech.com

[www.solarisbiotech.com](http://www.solarisbiotech.com)