IYGON

ND 100-80 Medical Tubing

Medical Products

From the ND Experts: The next generation in non-DEHP tubing

Features/Benefits

- High rigidity for kink resistance
- Meets USP Class VI criteria
- Compatible with blood
- Clear, with a glass-smooth inner bore
- Reduced flexibility for suction applications
- REACH compliant

Typical Applications

- Catheters
- Intravenous or intra-arterial infusion

Plasticizer

• TOTM

TYGON[®] ND 100-80 Tubing

TYGON[®] ND 100-80, one of the medical industry's first non-DEHP tubing offerings, has been specifically formulated in a higher durometer that makes it possible to produce smaller diameter tubings with a high resistance to kinking and collapse.

Characteristics

TYGON® ND 100-80 is available in five extremely small sizes for use in catheters, intravenous or intra-arterial infusion, and other surgical and hospital applications. It is rigid enough for easy handling, yet soft enough to reduce problems with puncturing.

Since TYGON® ND 100-80 is thermoplastic in nature, it can be fabricated into cannulas or protective sheath products using thermoforming and flaring techniques.

Based upon a unique chemistry developed by Saint-Gobain, TYGON® ND Series tubing was formulated to be Registration, Evaluation, Authorization and restriction of CHemical substances (REACH) compliant for DEHP. To ensure compliance, every individual compound lot is tested to ensure DEHP levels are <1000 ppm. During the extrusion process, individual product dimensions are maintained and monitored through in-line micrometers and off-line verification with computerized imaging equipment.

Consistent with many medical tubing requirements, ND 100-80 material can be effectively bonded/ welded using the following methods: heat, electronic (RF)/ultrasonic, solvent and adhesive. Factors to be considered when selecting the components include: security of the bond required, effect on the integrity of the materials to be joined, and presence of residues or extractables that may affect biocompatibility. When bonding procedures are not used, mechanical clamps are recommended to provide secure attachment.

Biocompatibility

TYGON® ND 100-80 is made from a biocompatible non-DEHP polymer material developed specifically for medical device needs. The product is USP Class VI approved and REACH compliant.



PERFORMANCE PLASTICS

Tel.: +52 (55) 6273-4046 / 6273-4048 / 2465-1016

TYGON[®] ND 100-80 Inventoried Sizes

Part Number	Inside Diameter		Outside Diameter		Wall Thickness		Coil Length	French Sizes	Needle Gauge
	Inches	mm	Inches	mm	Inches	mm	Feet		
AAD04091	.010	.25	.030	.76	.010	.25	500	3	30
AADo4103	.020	.51	.060	1.52	.020	.51	500	X 5 X	23
AADo4119	.030	.76	.090	2.29	.030	.76	500	7	20 – 21
AAD04127	.040	1.02	.070	1.78	.015	.38	500	6	18 – 20
AAD04133	.050	1.27	.090	2.29	.020	.51	500	10	17 – 18

TYGON[®] ND 100-80 Typical **Physical Properties**

Property	ASTM Method	Value or Rating	
Durometer Hardness, Shore A, 15 sec.	D2240	80	
Color		Clear	
Tensile Strength, psi (MPa)	D412	2,625 (18.1)	
Ultimate Elongation, %	D412	350	
Tear Resistance, lb-f/inch (kN/m)	D1004	275 (48.2)	
Specific Gravity	D792	1.21	
Water Absorption, % 24 hrs. @ 73°F (23°C)	D570	0.12	
Compression Set, Constant Deflection, % @ 158°F (70°C) for 22 hrs.	D395 Method B	59 8	
Maximum Recommended Operating Temp.,°F (°C)	8 83	185 (85)	
Brittleness by Impact Temp., °F (°C)	D746	-31 (-35)	
Tensile Modulus, @ 100% Elongation, psi (MPa)	D412	1,450 (10.0)	

Unless otherwise noted, all tests were conducted at room temperature (73°F). Values shown were determined on 0.075" thick extruded strip or 0.075" thick molded ASTM plaques or molded ASTM durometer buttons.

TYGON[®] ND 100-80 Characteristics

The biocompatibility of TYGON® ND 100-80 has been tested and found to be non-toxic in the following test protocols:

Opacity Color Color	Transparent		
FDA Approved for Food Contact	No		
NSF 51 Standard	No		
USP Class VI	Yes		

Sterilization Methods

Autoclavable	Yes		
Gas Co Co Co	Yes		
Radiation	Yes X		

TYGON[®] tubing is not intended for use as an implant material

Saint-Gobain Non-DEHP tubing contains DEHP levels of <1000 ppm.

TYGON[®] is a registered trademark.

SAINT-GOBAIN **PERFORMANCE PLASTICS**

Saint-Gobain Performance Plastics 2664 Gilchrist Road

Akron, OH 44305 Tel: 800-798-1554 Tel: (330) 798-9240 Fax: (330) 798-6968 **IMPORTANT:** It is the user's responsibility to ensure the suitability and safety of Saint-Gobain Performance Plastics tubing for all intended uses. Laboratory and clinical tests must be conducted in accordance with applicable regulatory requirements in order to determine the safety and effectiveness for use of tubing in any particular application

riod of 6 months from the date of first sale, Saint-Gobain Performance Corporation warrants this product to be free from defects in materials and proving defective nshin Our only obligation damage. direct or consequential, arising nisuse. or inability to use, this product, THIS WARRANTY IS IN LIEU INESS FOR PARTICULAR PURPOSE AND ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED. No deviation is authorized

Saint-Gobain Performance Plastics Corporation assumes no obligations or liability for any advice furnished by it, or for results obtained with respect to those products. All such advice is given and accepted at the buyer's risk

